LuxmedIQ

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LuxmedIQ

CE 0197

One-Step Drug of Abuse Test (Strip) Amphetamine (AMP 1000)

For self-testing and in vitro diagnostic use only

The LuxmedIQ One-Step Drug of Abuse Test is a rapid test for the qualitative detection of D-Amphetamine in human urine at a cutoff concentration of 1000 ng/mL. The test may yield preliminary positive results when prescription drugs are ingested at prescribed doses. It is not intended to distinguish between prescription use and abuse of any drug. There are no uniformly recognized cutoff concentration levels for any drug in urine. The test provides only preliminary test results. A more specific alternative chemical method must be used to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be exercised with any drug of abuse test result, particularly when the preliminary result is positive. For self-testing and in vitro diagnostic use only.

WHAT IS A ONE-STEP DRUG OF ABUSE TEST?


The LuxmedIQ One-Step Drug of Abuse Test is a rapid test for the qualitative detection of D-Amphetamine in human urine. The One-Step Drug of Abuse Test yields a positive result when then concentration of a drug and/or its metabolite in urine is at or exceeds its cutoff concentration.

WHAT IS THE CUTOFF VALUE AND APPROXIMATE DETECTION TIME?

Drug (Identifier)	Cutoff Level	Minimum Detection Time	Maximum Detection Time
D-Amphetamine	1,000 ng/mL	4-6 hours	2-3 days

PRINCIPLE

The LuxmedIQ One-Step Drug of Abuse Test is an immunoassay. During testing, the urine specimen migrates upward on the test strip. A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip, while a drug-negative urine specimen will generate a line in the test line region. A colored line will always appear in the control line region, indicating that the proper specimen volume has been added.

REAGENTS

The test contains a membrane strip coated with drug-protein conjugates (purified bovine albumin) on the test line, a goat polyclonal antibody against gold-protein conjugate at the control line, and a dye pad which contains colloidal gold particles coated with mouse monoclonal antibody specific to cocaine.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use.
- For external use only.
- For single use. Discard after first use.
- Do not use the test if the pouch is punctured or not well sealed.
- Do not use after the expiration date.
- Keep out of the reach of children.
- The used urine sample cup and the urine specimen should be discarded according to federal, state and local regulations.
- The used test device and the urine specimen should be discarded according to federal, state and local regulations.

CONTENT OF THE PACKAGE

INCLUDED IN THE PACKAGE:

- User instruction
- Test strip (inside a foil pouch)
- Desiccants

NOT INCLUDED IN THE PACKAGE:

- Watch, timer or clock
- Disposable gloves
- Collection cup for the urine sample

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 39°F - 86°F (4°C - 30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. If the test remains unused for more than 10 minutes after removal from the sealed pouch, results may be erroneous. Keep away from direct sunlight, moisture, and heat. **DO NOT FREEZE.** Do not use beyond the expiration date.

WHEN TO COLLECT THE URINE SAMPLE FOR THE TEST?

Urine from any time of the day can be used.


HOW TO COLLECT THE URINE SAMPLE?

- When you are ready to test, remove the test strip from the sealed foil pouch.
- Notice the colored tape on the strip indicates the name of the drug you are testing for.
- Fill a collection cup with a fresh urine sample. **DO NOT FREEZE.**

HOW TO DO THE TEST?

- Remove the test strip from the sealed pouch and use it as soon as possible.
- Immerse the absorbent tip of the test strip into the urine for a minimum of 15 seconds. Do not immerse the test strip beyond the MAX (maximum) line. Take the test strip out of the urine sample and lay it flat on a clean non-absorbent surface.
- Wait for 5 minutes (start timing immediately after the strip is taken out of the urine sample) and read the results. **DO NOT** read the results after more than 5 minutes.

NOTE: Results read after over 5 minutes may not be accurate and should not be used.

READING THE RESULTS

PRELIMINARY POSITIVE (+)

If a line appears in the C - Control area but NO line appears in the T - Drug Test area, then this indicates a preliminary positive result for the corresponding drug.

NEGATIVE (-)

If a line appears in both the C - Control and T - Drug Test area, then this indicates a negative result for the corresponding drug, regardless of how dark or light the line may appear.

INVALID

If at 5 minutes, NO line appears in the C - Control area, then the result is invalid.

NOTE: Each test strip must be looked at individually. Each line may vary in color and darkness. (**DO NOT** compare lines within the same test strip or between different test strips).

A positive test result does not always mean a person took illegal drugs and a negative test result does not always mean a person did not take illegal drugs. There are several factors that influence the reliability of drug tests. Certain drugs of abuse tests are more accurate than others.

IMPORTANT: The result you obtained is called preliminary for a reason. The sample must be tested by laboratory to determine if a drug of abuse is actually present.

WHAT IS A FALSE POSITIVE TEST?

The definition of a false positive test would be an instance where the test result from the LuxmedIQ One-Step Drug of Abuse Test is positive, even though the initial target drug is not present in the sample. The most common causes of a false positive test are cross reactants. Certain foods and medicines, dietary plan drugs and nutritional supplements may also cause a false positive test result with this product.

WHAT IS A FALSE NEGATIVE TEST?

The definition of a false negative test is that the initial target drug is present in the urine sample but is not detected by the LuxmediQ One-Step Drug of Abuse Test. A diluted sample or a sample contaminated with a substance could cause a false negative result.

TEST LIMITATIONS

- 1. The LuxmediQ One-Step Drug of Abuse Test provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
- 2. There is the possibility that interfering substances in the urine specimen may cause erroneous results.
- 3. Substances, such as bleach and/or alum, in urine samples may produce erroneous results.
- 4. A positive result does not indicate intoxication, the concentration of the drug in the urine, or the route of drug administration.
- 5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when the drug is present but below the cutoff level of the test.
- 6. The test does not distinguish between drugs of abuse and certain medications.
- 7. A positive test result may be obtained from certain foods or food supplements.

QUALITY CONTROL

If you work in a laboratory, you should perform quality control testing and you should read this section. A procedural control is included in the test. A color line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance. Quality control testing should be done with each new lot and each new shipment. It should be done every thirty days to check storage. Please contact our Technical Support at service@luxmediq.de for controls that work with the test.

PERFORMANCE CHARACTERISTICS

80 clinical urine specimens were analyzed by GC/MS and by the LuxmediQ One-Step Drug of Abuse Test. Each test was read by three viewers. Samples were divided by concentration into five categories: drug-free, less than half the cutoff, near cutoff negative, near cutoff positive, and high positive. Results were as follows:

ACCURACY - D-AMPHETAMINE

Viewer A:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	13	26
Negative	10	10	19	1	0

% agreement among positives is 97.5%
% agreement among negatives is 97.5%

Viewer B:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	12	26
Negative	10	10	19	2	0

% agreement among positives is 95%
% agreement among negatives is 97.5%

Viewer C:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	0	13	26
Negative	10	10	20	1	0

% agreement among positives is 97.5%
% agreement among negatives is 100%

From the results of the above tables, the total results are shown as below for D-Amphetamine:
The average positive agreement is 96.7%.
The average negative agreement is 98.3%.

PRECISION AND SENSITIVITY - D-AMPHETAMINE

Lot 1

Approximate Concentration of Sample (ng/mL)	Number of Determinations	Results Negative/Positive
0	50	50/0
250	50	50/0
500	50	50/0
750	50	50/0
1,000	50	2/48
1,250	50	0/50
1,500	50	0/50
1,750	50	0/50
2,000	50	0/50

Lot 2

Approximate Concentration of Sample (ng/mL)	Number of Determinations	Results Negative/Positive
0	50	50/0
250	50	50/0
500	50	50/0
750	50	50/0
1,000	50	3/47
1,250	50	0/50
1,500	50	0/50
1,750	50	0/50
2,000	50	0/50

Lot 3

Approximate Concentration of Sample (ng/mL)	Number of Determinations	Results Negative/Positive
0	50	50/0
250	50	50/0
500	50	50/0
750	50	50/0
1,000	50	1/49
1,250	50	0/50
1,500	50	0/50
1,750	50	0/50
2,000	50	0/50

SPECIFICITY AND CROSS REACTIVITY

To test the specificity of the test, the test device was used to test D-Amphetamine, drug metabolites and other components of the same class that are likely to be present in urine. All the components were added to drug-free normal human urine. The following structurally related compounds produced positive results with the test when tested at levels equal to or greater than the concentrations listed below.

D-Amphetamine	Result
(D-Amphetamine, Cutoff = 1000ng/mL)	Positive at 1,000ng/mL
L-Amphetamine	Positive at 100,000ng/mL
DL-Amphetamine	Positive at 500ng/mL
(+/-)-3,4-Methylenedioxymphetamine (MDA)	Positive at 1,300ng/mL
Phentermine	Positive at 100,000ng/mL
Apomorphine	Positive at 50,000ng/mL
β-Phenylethylamine	Positive at 25,000ng/mL
Tyramine	Positive at 10,000ng/mL
Tryptamine	Positive at 25,000ng/mL
D-Methamphetamine	> 100,000
L-Methamphetamine	> 100,000
Ephedrine	> 100,000

3,4-Methylenedioxyethylamphetamine (MDE)	> 100,000	(1R,2S)-(-)-Ephedrine	Oxolinic Acid	Zomepirac
--	-----------	-----------------------	---------------	-----------

EFFECT OF URINARY SPECIFIC GRAVITY

Urine samples of normal, high, and low specific gravity ranges (1.000-1.035) were spiked with drugs at 25% below and 25% above cutoff levels respectively. The Lux-medIQ One-Step Drug of Abuse Test was tested using twelve drug-free urine and spiked urine samples. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

EFFECT OF URINARY PH

The pH of an aliquot of a negative urine pool was adjusted to pH ranges of 4.0 - 9.0 and spiked with drugs at 25% below and 25% above cut-off levels. The spiked, pH-adjusted urine was tested with the LuxmedIQ One-Step Drug of Abuse Test. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

NON CROSS-REACTING COMPOUNDS - D-AMPHETAMINE

4-Acetamidophenol	L-Ephedrine	Oxycodone
Acetophenetidin	(-)-ψ-Ephedrine	Oxymetazoline
N-Acetylprocainamide	Erythromycin	Papaverine
Acetylsalicylic Acid	β-Estradiol	Penicillin-G
Aminopyrine	Estrone-3-Sulfate	Pentazocine
Amitriptyline	Ethyl-p-Aminobenzoate	Pentobarbital
Amobarbital	Fenfluramine	Perphenazine
Amoxicillin	Fenoprofen	Phencyclidine
Ampicillin	Furosemide	Phenelzine
Ascorbic Acid	Gentisic Acid	Phenobarbital
Aspartame	Hemoglobin	Phenytoin
Atropine	Hydralazine	L-Phenylephrine
Benzilic Acid	Hydrochlorothiazide	Phenylpropanolamine
Benzoic Acid	Hydrocodone	Prednisolone
Benzoylcegonine	Hydrocortisone	Prednisone
Bilirubin	o-Hydroxyhippuric Acid	Procaine
Brompheniramine	3-Hydroxytyramine	Promazine
Caffeine	Ibuprofen	Promethazine
Cannabidiol	Imipramine	DL-Propranolol
Cannabinalol	(-)-Isoproterenol	D-Propoxyphene
Chloral Hydrate	Isoxsuprine	Quinidine
Chloramphenicol	Ketamine	Quinine
Chlordiazepoxide	Ketoprofen	Ranitidine
Chlorothiazide	Labetalol	Salicylic Acid
(±)-Chlorpheniramine	Levorphanol	Secobarbital
Chlorpromazine	Loperamide	Sulfamethazine
Chloroquine	Maprotiline	Sulindac
Cholesterol	Meperidine	Temazepam
Clomipramine	Meprobamate	Tetracycline
Clonidine	Methadone	Tetrahydrocortisone
Cocaine Hydrochloride	Methylphenidate	Tetrahydrozoline
Codeine	Morphine 3-β-D-Glu-curonide	Δ ⁸ -THC-COOH
Cortisone	Nalidixic Acid	Thebaine
(-)-Cotinine	Naloxone	Thiamine
Creatinine	Naltrexone	Thioridazine
Deoxycorticosterone	Naproxen	DL-Thyroxine
Dextromethorphan	Niacinamide	Tolbutamide
Diazepam	Nifedipine	Triamterene
Diclofenac	Norcodeine	Trifluoperazine
Diflunisal	Norethindrone	Trimethoprim
Digoxin	D-Norpropoxyphene	Trimipramine
Diphenhydramine	Noscapine	Tryptamine
Doxylamine	DL-Octopamine	DL-Tyrosine
Ecgonine Hydrochloride	Oxalic Acid	Uric Acid
Ecgonine Methyl Ester	Oxazepam	Verapamil

LAY USER

A lay user study was performed at three intended user sites with 140 laypersons for each drug device. They had diverse educational and professional backgrounds and ranged in age from 21 to >50. Urine samples were prepared at the following concentrations: negative, +/-75%, +/-50%, +/-25% of the cutoff by spiking drug(s) into drug free-pooled urine specimens. The concentrations of the samples were confirmed by GC/MS. Each sample was aliquoted into individual containers and blind-labeled. Each participant was provided with the package insert, 1 blind labeled sample and a test strip. The results are summarized below.

% of Cutoff	Number of Samples	D-Amphet-amine Con- centration by GC/MS (ng/mL)	No. of Positive	No. of Negative	The Percentage Agreement (%)
-100%Cutoff	20	0	0	20	100%
-75%Cutoff	20	246	0	20	100%
-50%Cutoff	20	492	0	20	100%
-25%Cutoff	20	738	2	18	90%
+25%Cutoff	20	1,268	19	1	95%
+50%Cutoff	20	1,521	20	0	100%
+75%Cutoff	20	1,775	20	0	100%

BIBLIOGRAPHY OF SUGGESTED READINGS

- Stewart DJ, Inaba T, Lucassen M, Kalow W. Clin. Pharmacol. Ther. April 1979; 25 ed: 464, 264-8.
- Ambre J. J. Anal. Toxicol. 1985; 9:241.
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

SYMBOLS

	Do not reuse		In vitro diagnostic medical device
	Temperature limitations		Consult instructions before use
	Manufacturer		Authorized Representative in the European Union
	Batch code		Use by
	CE Mark		Catalog number
	Contains sufficient for <n> tests		

DISTRIBUTED BY:

LuxmediQ Deutschland UG
Altonaer Straße 27, 10555 Berlin, Germany

W.H.P.M. Bioresearch & Technology Co., Ltd.
No. 2 Zhongxin Street, Louzizhuang, Jinzhan Xiang, Chaoyang District, Beijing,
10018 P.R. China
Tel: (+86)-10-84391888
Fax: (+86)-10-84391888-117
www.whpm.com.cn

Medical Device Safety Service GmbH
Schiffgraben 41, 30175 Hannover, Germany
Tel: (+49)-511-6262-8630
Fax: (+49)-511-6262 8633
www.mdss.com

Effective date: 2017-09-01
Y0311153301

FRAGEN | QUESTIONS?

E-mail: service@luxmediq.de

Tel: +49 6401 9699 390

Web: www.luxmediq.de