LuxmedIQ

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LuxmedIQ

CE 0197

One-Step Drug of Abuse Test (Strip) Ecstasy (MDMA 500)

For self-testing and in vitro diagnostic use only

The LuxmedIQ One-Step Drug of Abuse Test is a rapid test for the qualitative detection of Methylenedioxyamphetamine in human urine at a cutoff concentration of 150 ng/mL.

The test may yield preliminary positive results when prescription drugs are ingested at prescribed doses. It is not intended to distinguish between prescription use and abuse of any drug. There are no uniformly recognized cutoff concentration levels for any drug in urine. The test provides only preliminary test results. A more specific alternative chemical method must be used to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be exercised with any drug of abuse test result, particularly when the preliminary result is positive.

For self-testing and in vitro diagnostic use only.

WHAT IS A ONE-STEP DRUG OF ABUSE TEST?


The LuxmedIQ One-Step Drug of Abuse Test is a rapid test for the qualitative detection of Methylenedioxyamphetamine in human urine. The One-Step Drug of Abuse Test yields a positive result when then concentration of a drug and/or its metabolite in urine is at or exceeds its cutoff concentration.

WHAT IS THE CUTOFF VALUE AND APPROXIMATE DETECTION TIME?

Drug (Identifier)	Cutoff Level	Minimum Detection Time	Maximum Detection Time
Methylenedioxy-methamphetamine	500 ng/mL	2-7 hours	2-4 days

PRINCIPLE

The LuxmedIQ One-Step Drug of Abuse Test is an immunoassay. During testing, the urine specimen migrates upward on the test strip. A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip, while a

drug-negative urine specimen will generate a line in the test line region. A colored line will always appear in the control line region, indicating that the proper specimen volume has been added.

REAGENTS

The test contains a membrane strip coated with drug-protein conjugates (purified bovine albumin) on the test line, a goat polyclonal antibody against gold-protein conjugate at the control line, and a dye pad which contains colloidal gold particles coated with mouse monoclonal antibody specific to cocaine.

WARNINGS AND PRECAUTIONS

- For in vitro diagnostic use.
- For external use only.
- For single use. Discard after first use.
- Do not use the test if the pouch is punctured or not well sealed.
- Do not use after the expiration date.

- Keep out of the reach of children.
- The used urine sample cup and the urine specimen should be discarded according to federal, state and local regulations.
- The used test device and the urine specimen should be discarded according to federal, state and local regulations.

CONTENT OF THE PACKAGE

INCLUDED IN THE PACKAGE:

- User instruction
- Test strip (inside a foil pouch)
- Desiccants

NOT INCLUDED IN THE PACKAGE:

- Watch, timer or clock
- Disposable gloves
- Collection cup for the urine sample

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 39°F - 86°F (4°C - 30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. If the test remains unused for more than 10 minutes after removal from the sealed pouch, results may be erroneous. Keep away from direct sunlight, moisture, and heat. **DO NOT FREEZE.** Do not use beyond the expiration date.

WHEN TO COLLECT THE URINE SAMPLE FOR THE TEST?

Urine from any time of the day can be used.

HOW TO COLLECT THE URINE SAMPLE?

- When you are ready to test, remove the test strip from the sealed foil pouch.
- Notice the colored tape on the strip indicates the name of the drug you are testing for.
- Fill a collection cup with a fresh urine sample. **DO NOT FREEZE.**

HOW TO DO THE TEST?

- Remove the test strip from the sealed pouch and use it as soon as possible.
- Immerse the absorbent tip of the test strip into the urine for a minimum of 15 seconds. Do not immerse the test strip beyond the MAX (maximum) line. Take the test strip out of the urine sample and lay it flat on a clean non-absorbent surface.
- Wait for 5 minutes (start timing immediately after the strip is taken out of the urine sample) and read the results. **DO NOT** read the results after more than 5 minutes.

NOTE: Results read after over 5 minutes may not be accurate and should not be used.

READING THE RESULTS

PRELIMINARY POSITIVE (+)

If a line appears in the C - Control area but NO line appears in the T - Drug Test area, then this indicates a preliminary positive result for the corresponding drug.

NEGATIVE (-)

If a line appears in both the C - Control and T - Drug Test area, then this indicates a negative result for the corresponding drug, regardless of how dark or light the line may appear.

INVALID

If at 5 minutes, NO line appears in the C - Control area, then the result is invalid.

NOTE: Each test strip must be looked at individually. Each line may vary in color and darkness. **(DO NOT)** compare lines within the same test strip or between different test strips.

A positive test result does not always mean a person took illegal drugs and a negative test result does not always mean a person did not take illegal drugs. There are several factors that influence the reliability of drug tests. Certain drugs of abuse tests are more accurate than others.

IMPORTANT: The result you obtained is called preliminary for a reason. The sample must be tested by laboratory to determine if a drug of abuse is actually present.

WHAT IS A FALSE POSITIVE TEST?

The definition of a false positive test would be an instance where the test result from the LuxmedIQ One-Step Drug of Abuse Test is positive, even though the initial target drug is not present in the sample. The most common causes of a false positive test are cross reactants. Certain foods and medicines, dietary plan drugs and nutritional supplements may also cause a false positive test result with this product.

WHAT IS A FALSE NEGATIVE TEST?

The definition of a false negative test is that the initial target drug is present in the urine sample but is not detected by the LuxmedIQ One-Step Drug of Abuse Test. A diluted sample or a sample contaminated with a substance could cause a false negative result.

TEST LIMITATIONS

1. The LuxmedIQ One-Step Drug of Abuse Test provides only a qualitative, preliminary analytical result. A secondary analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method.
2. There is the possibility that interfering substances in the urine specimen may cause erroneous results.
3. Substances, such as bleach and/or alum, in urine samples may produce erroneous results.
4. A positive result does not indicate intoxication, the concentration of the drug in the urine, or the route of drug administration.
5. A negative result may not necessarily indicate drug-free urine. Negative results can be obtained when the drug is present but below the cutoff level of the test.
6. The test does not distinguish between drugs of abuse and certain medications.
7. A positive test result may be obtained from certain foods or food supplements.

QUALITY CONTROL

If you work in a laboratory, you should perform quality control testing and you should read this section. A procedural control is included in the test. A color line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique. Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and to verify proper test performance. Quality control testing should be done with each new lot and each new shipment. It should be done every thirty days to check storage. Please contact our Technical Support at service@luxmediq.de for controls that work with the test.

PERFORMANCE CHARACTERISTICS

80 clinical urine specimens were analyzed by GC/MS and by the LuxmedIQ One-Step Drug of Abuse Test. Each test was read by three viewers. Samples were divided by concentration into five categories: drug-free, less than half the cutoff, near cutoff negative, near cutoff positive, and high positive. Results were as follows:

ACCURACY - METHYLENEDIOXYMETHAMPHETAMINE

Viewer A:

WHPM Result	Drug-Free	Less Than Half The Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	14	26
Negative	10	10	19	0	0

% agreement among positives is 100%
% agreement among negatives is 97.5%

Viewer B:

WHPM Result	Drug-Free	Less Than Half The Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	0	13	26
Negative	10	10	20	1	0

% agreement among positives is 97.5%
% agreement among negatives is 100%

Viewer C:

WHPM Result	Drug-Free	Less Than Half The Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	14	26
Negative	10	10	19	0	0

% agreement among positives is 100%
% agreement among negatives is 97.5%

From the results of the above tables, the total results are shown as below for Methylenedioxyamphetamine:
The average positive agreement is 99.2%.
The average negative agreement is 98.3%.

PRECISION AND SENSITIVITY - METHYLENEDIOXYMETHAMPHETAMINE

Lot 1

Approximate Concentration of Sample (ng/mL)	Number of Determinations	Results Negative/Positive
0	50	50/0
125	50	50/0
250	50	50/0
375	50	50/0
500	50	3/47
625	50	0/50
750	50	0/50
875	50	0/50
1,000	50	0/50

Lot 2

Approximate Concentration of Sample (ng/mL)	Number of Determinations	Results Negative/Positive
0	50	50/0
125	50	50/0
250	50	50/0
375	50	50/0
500	50	2/48
625	50	0/50
750	50	0/50
875	50	0/50
1,000	50	0/50

Lot 3

Approximate Concentration of Sample (ng/mL)	Number of Determinations	Results Negative/Positive
0	50	50/0
125	50	50/0
250	50	50/0
375	50	50/0
500	50	2/48
625	50	0/50
750	50	0/50
875	50	0/50
1,000	50	0/50

SPECIFICITY AND CROSS REACTIVITY

To test the specificity of the test, the test device was used to test Methylenedioxyamphetamine, drug metabolites and other components of the same class that are likely to be present in urine. All the components were added to drug-free normal human urine. The following structurally related compounds produced positive results with the test when tested at levels equal to or greater than the concentrations listed below.

Methylenedioxymethamphetamine	Result
(Methylenedioxy-methamphetamine, Cutoff = 500ng/mL)	Positive at 500ng/mL
3,4-Methylenedioxyamphetamine HCl (MDA)	Positive at 8,000ng/mL
3,4-Methylenedioxyethylamphetamine (MDEA)	Positive at 1,000ng/mL
(-)- ψ -Ephedrine	Positive at 40,000ng/mL
D-Methamphetamine	Negative at $\leq 10^3$ ng/mL
D-Amphetamine	Negative at $\leq 10^3$ ng/mL
L-Amphetamine	Negative at $\leq 10^3$ ng/mL
L-Methamphetamine	Negative at $\leq 10^3$ ng/mL

EFFECT OF URINARY SPECIFIC GRAVITY

Urine samples of normal, high, and low specific gravity ranges (1.000-1.035) were spiked with drugs at 25% below and 25% above cutoff levels respectively. The LuxmedQ One-Step Drug of Abuse Test was tested using twelve drug-free urine and spiked urine samples. The results demonstrate that varying ranges of urinary specific gravity do not affect the test results.

EFFECT OF URINARY pH

The pH of an aliquot of a negative urine pool was adjusted to pH ranges of 4.0 - 9.0 and spiked with drugs at 25% below and 25% above cut-off levels. The spiked, pH-adjusted urine was tested with the LuxmedQ One-Step Drug of Abuse Test. The results demonstrate that varying ranges of pH do not interfere with the performance of the test.

NON CROSS-REACTING COMPOUNDS - MDMA

4-Acetamidophenol	L-Epinephrine	Pentobarbital
Acetophenetidin	Erythromycin	Perphenazine
N-Acetylprocainamide	β -Estradiol	Phencyclidine
Acetylsalicylic Acid	Estrone-3-Sulfate	Phenelzine
Aminopyrine	Ethyl-p-Aminobenzoate	Phenobarbital
Amitriptyline	Fenoprofen	Phentermine
Amobarbital	Furosemide	Trans-2-Phenylcyclopropylamine Hydrochloride
Amoxicillin	Gentisic Acid	L-Phenylephrine
Ampicillin	Hemoglobin	β -Phenylethylamine
L-Ascorbic Acid	Hydralazine	Phenylpropanolamine
Apomorphine	Hydrochlorothiazide	Prednisolone
Aspartame	Hydrocodone	Prednisone
Atropine	Hydrocortisone	Procaine
Benzilic Acid	o-Hydroxyhippuric Acid	Promazine
Benzoic Acid	3-Hydroxytyramine	Promethazine
Benzoylcegonine	Ibuprofen	DL-Propranolol
Bilirubin	Imipramine	D-Propoxyphene
(\pm)-Brompheniramine	Iproniazid	D-Pseudoephedrine
Bupropion	(\pm)-Isoproterenol	Quinacrine
Caffeine	Isoxsuprine	Quinidine
Cannabidiol	Ketamine	Quinine
Cannabinol	Ketoprofen	Ranitidine
Chloral Hydrate	Labetalol	Salicylic Acid
Chloramphenicol	Levorphanol	Secobarbital
Chlordiazepoxide	Loperamide	Serotonin (5-Hydroxy-tyramine)
Chlorothiazide	Maprotiline	Sulfamethazine
(\pm)-Chlorpheniramine	Meperidine	Sulindac
Chlorpromazine	Meprobamate	Sustiva
Chloroquine Methylphenidate	Methodone	Temazepam
Cholesterol	Morphine-3- β -D-Glucuronide	Tetracycline
Clomipramine	Morphine Sulfate	Tetrahydrocortisone 3- β -D-Glucuronide)
Clonidine	Nalidixic Acid	Tetrahydrozoline
Cocaethylene	Naloxone	Thebaine

Cocaine Hydrochloride	Naltrexone	Theophylline
Codeine	Naproxen	Thiamine
Cortisone	Niacinamide	Thioridazine
(-)-Cotinine	Nifedipine	Tolbutamide
Creatinine	Nimesulide	Trazodone
Deoxycorticosterone	Norcodeine	Triamterene
Dextromethorphan	Norethindrone	DL-Tyrosine
Diclofenac	D-Norpropoxyphene	Trifluoperazine
Diazepam	Noscapine	Trimethoprim
Diffunisal	DL-Octopamine	Trimipramine
Digoxin	Oxalic Acid	Tryptamine
Dicyclomine	Oxazepam	DL-Tryptophan
Diphenhydramine	Oxolinic Acid	Tyramine
5,5-Diphenylhydantoin	Oxycodone	Uric Acid
Doxylamine	Oxymetazoline	Verapamil
Ecgonine Hydrochloride	Papaverine	Zomepirac
Ecgonine Methyl Ester	Penicillin-G	(1R,2S)-(-)-Ephedrine
Pentazocine Hydrochloride		

LAY USER


A lay user study was performed at three intended user sites with 140 laypersons for each drug device. They had diverse educational and professional backgrounds and ranged in age from 21 to >50. Urine samples were prepared at the following concentrations: negative, +/-75%, +/-50%, +/-25% of the cutoff by spiking drug(s) into drug free-pooled urine specimens. The concentrations of the samples were confirmed by GC/MS. Each sample was aliquoted into individual containers and blind-labeled. Each participant was provided with the package insert, 1 blind labeled sample and a test strip. The results are summarized below.

% of Cutoff	Number of Samples	MDMA Concentration by GC/MS (ng/mL)	No. of Positive	No. of Positive	The Percentage Agreement (%)
-100% Cutoff	20	0	0	20	100%
-75% Cutoff	20	115	0	20	100%
-50% Cutoff	20	237	0	20	100%
-25% Cutoff	20	358	1	19	95%
+25% Cutoff	20	598	20	0	100%
+50% Cutoff	20	755	20	0	100%
+75% Cutoff	20	912	20	0	100%

BIBLIOGRAPHY OF SUGGESTED READINGS

- Stewart DJ, Inaba T, Lucassen M, Kalow W. Clin. Pharmacol. Ther. April 1979; 25 ed: 464, 264-8.
- Ambre J. J. Anal. Toxicol. 1985; 9:241.
- Hawks RL, CN Chiang. Urine Testing for Drugs of Abuse. National Institute for Drug Abuse (NIDA), Research Monograph 73, 1986.

SYMBOLS

	Do not reuse		In vitro diagnostic medical device
	Temperature limitations		Consult instructions before use
	Manufacturer		Authorized Representative in the European Union
	Batch code		Use by
	CE Mark		Catalog number
	Contains sufficient for <n> tests		

DISTRIBUTED BY:

LuxmedIQ Deutschland UG
Altonaer Straße 27, 10555 Berlin, Germany

W.H.P.M. Bioresearch & Technology Co., Ltd.
No. 2 Zhongxin Street, Louzizhuang, Jinzhan Xiang, Chaoyang District, Beijing,
10018 P.R. China
Tel: (+86)-10-84391888
Fax: (+86)-10-84391888-117
www.whpm.com.cn

Medical Device Safety Service GmbH
Schiffgraben 41, 30175 Hannover, Germany
Tel: (+49)-511-6262-8630
Fax: (+49)-511-6262-8633
www.mdss.com

Effective date: 2017-09-01
Y0311153301

FRAGEN | QUESTIONS?

E-mail: service@luxmediq.de

Tel: +49 6401 9699 390

WhatsApp: +49 6401 9699 390

Web: www.luxmediq.de