LuxmedIQ

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## BIBLIOGRAPHIE

1. Moolchan, E. et al, "Saliva and Plasma Testing for Drugs of Abuse: Comparison of the Disposition and Pharmacological Effects of Cocaine", Addiction Research Center, IRP, NIDA, NIH, Baltimore, MD. As presented at the SOFT-TIAFT meeting October 1998.
2. Kim, I, et al, "Plasma and oral fluid pharmacokinetics and pharmacodynamics after oral codeine administration", Clin Chem, 2002 Sept; 48 (9), pp 1486-96.
3. Schramm, W. et al, "Drugs of Abuse in Saliva: A Review," J Anal Tox, 1992 Jan-Feb; 16 (1), pp 1-9.
4. McCarron, MM, et al, "Detection of Phencyclidine Usage by Radioimmunoassay of Saliva," J Anal Tox. 1984 Sep-Oct; 8 (5), pp 197-201.
5. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735.

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LuxmedIQ

# LuxmedIQ One-Step Drug of Abuse Test (Oral Fluid)

## Multi 11 + ALC

For self-testing use only

The LuxmedIQ One-Step Drug of Abuse Test Multi 11 + ALC is a rapid test for the qualitative detection of Amphetamine, Barbiturates, Benzodiazepines, Cocaine, Methamphetamines, Methylenedioxy-methamphetamine, Methadone, Opiates, Oxycodone, Phencyclidine, Marijuana, and Alcohol in human oral fluid.

**NOTE:** Please read these instructions carefully, because they contain important information on what should be considered when performing the test. Familiarize yourself with these instructions before using the contents of the package.

### INTENDED USE

#### For self-testing use only.

The Oral Fluid Drug Screen Device is a lateral flow chromatographic immunoassay for the qualitative detection of multiple drugs or drug metabolites in oral fluids at the following cut-off concentrations:

Test (Identifier)	Calibrator	Cutoff
Amphetamine (AMP)	D-Amphetamine	50ng/mL
Barbiturates (BARB)	Butalbital	50ng/mL
Benzodiazepines (BZO)	Oxazepam	50ng/mL
Cocaine (COC)	Benzoylecgonine	20ng/mL
Methamphetamines (mAMP)	D-Methamphetamine	50ng/mL
Methylenedioxy-methamphetamine (MDMA)	(±)-3,4-Methylenedioxy-methamphetamine	50ng/mL
Methadone (MTD)	Methadone	30ng/mL
Opiates (OPI)	Morphine	40ng/mL
Oxycodone (OXY)	Oxycodone	50ng/mL
Phencyclidine (PCP)	Phencyclidine	10ng/mL
Marijuana (THC)	11-nor- $\Delta^9$ -THC-9-COOH	12ng/mL
Alcohol (ALCO)	Alcohol	>0.02% B.A.C.

This assay provides only a preliminary analytical test result. A more specific alternate chemical method must be used in order to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) and Liquid chromatography/tandem mass spectrometry (LC/MS/MS) are the preferred confirmatory methods. Professional judgment should be applied to any drug of abuse test result, particularly when preliminary results are indicated.

### SUMMARY AND EXPLANATION OF THE TEST

The Oral Fluid Drug Screen Device is a rapid, oral fluid screening test that can be performed without the use of an instrument. The test utilizes monoclonal antibodies to selectively detect elevated levels of specific drugs in human oral fluid.

#### AMPHETAMINE (AMP)

Amphetamine is a sympathomimetic amine with therapeutic indications. The drug is often self-administered by nasal inhalation or oral ingestion. Depending on the route of administration, Amphetamine can be detected in oral fluid as early as 5-10 minutes and up to 72 hours after use. The Amphetamine assay contained within the Oral Fluid Drug Screen Device yields a positive result when the Amphetamine concentration in oral fluid exceeds 50 ng/mL.

#### BARBITURATES (BARB)

Barbiturates are CNS depressants. They are used therapeutically as sedatives, hypnotics, and anticonvulsants. Barbiturates are almost always taken orally as capsules or tablets. The effects resemble those of intoxication with alcohol. Chronic use of barbiturates leads to tolerance and physical dependence. Short-acting barbiturates taken at 400 mg/day for 2-3 months can produce a clinically significant degree of physical dependence. Withdrawal symptoms experienced during periods of drug abstinence can be severe enough to cause death. Only a small amount (less than 5%) of most barbiturates are excreted unaltered in the urine. The approximate detection time limits

for barbiturates are:

Short acting (e.g. Secobarbital) 100 mg PO (oral) 4.5 days

Long acting (e.g. Phenobarbital) 400 mg PO (oral) 7 days

The Barbiturates assay contained within the Oral Fluid Drug Screen Device yields a positive result when the Butalbital concentration in oral fluid exceeds 50 ng/mL.

#### BENZODIAZEPINES (BZO)

Benzodiazepines are frequently prescribed sedative and hypnotic drug for the symptomatic treatment of anxiety, insomnia, sleep and seizure disorders. Most Benzodiazepines are extensively metabolized in the liver and excreted in the urine and saliva as metabolites. Chronic abuse may increase the risk of physical dependence and may result in intoxication, drowsiness, and muscle relaxation. Oxazepam is the major metabolic product of Benzodiazepines.

The Benzodiazepines assay contained within the Oral Fluid Drug Screen Device yields a positive result when the concentration of Oxazepam in oral fluids exceeds 50 ng/mL.

#### COCAINE (COC)

Cocaine is a potent central nervous system (CNS) stimulant and a local anesthetic derived from the coca plant (Erythroxylum coca). The drug is often self-administered by nasal inhalation, intravenous injection and free-base smoking. Depending on the route of administration, cocaine and metabolites benzoylecgonine and ecgonine methyl ester can be detected in oral fluid as early as 5-10 minutes following use. Cocaine and benzoylecgonine can be detected in oral fluids for up to 24 hours after use.

The Cocaine assay contained within the Oral Fluid Drug Screen Device yields a positive result when the cocaine metabolite in oral fluid exceeds 20 ng/mL.

#### METHAMPHETAMINE (MAMP)

Methamphetamine is a potent stimulant chemically related to amphetamine but with greater CNS stimulation properties. The drug is often self-administered by nasal inhalation, smoking or oral ingestion. Depending on the route of administration, methamphetamine can be detected in oral fluid as early as 5-10 minutes and up to 72 hours after use.

The Methamphetamine assay contained within the Oral Fluid Drug Screen Device yields a positive result when the Methamphetamine concentration in oral fluid exceeds 50 ng/mL.

#### METHYLENEDIOXYMETHAMPHETAMINE (MDMA)

Methylenedioxy-methamphetamine (ecstasy) is a designer drug first synthesized in 1914 by a German drug company for the treatment of obesity. Those who take the drug frequently report adverse effects, such as increased muscle tension and sweating. MDMA is not clearly a stimulant, although it has in common with amphetamine drugs, a capacity to increase blood pressure and heart rate. MDMA does produce some perceptual changes in the form of increased sensitivity to light, difficulty in focusing, and blurred vision in some users. Its mechanism of action is thought to be via release of the neurotransmitter serotonin. MDMA may also release dopamine, although the general opinion is that this is a secondary effect of the drug (Nichols and Oberlander, 1990). The most pervasive effect of MDMA, occurring in virtually all people who took a reasonable dose of the drug, was to produce a clenching of the jaws.

The Methylenedioxy-methamphetamine assay contained within the Oral Fluid Drug Screen Device yields a positive result when the Methylenedioxy-methamphetamine concentration in oral fluid exceeds 50 ng/mL.

#### METHADONE (MTD)

Methadone is a narcotic analgesic prescribed for the management of moderate to severe pain and for the treatment of opiate dependence (heroin, Vicodin, Percocet, morphine). The pharmacology of oral methadone is very different from IV methadone. Oral methadone is partially stored in the liver for later use. IV methadone acts more like heroin. In most states you must go to a pain clinic or a methadone maintenance clinic to be prescribed methadone. Methadone is a long-acting pain reliever producing relief that lasts from two to forty-eight hours. Ideal methadone frees the client from the pressures of obtaining illegal heroin, from the dangers of injection, and from the emotional roller coaster that most opiates produce. Methadone, if taken for long periods and at large doses, can lead to a very long withdrawal period. The withdrawals from methadone are more prolonged and troublesome than those provoked by heroin cessation, yet the substitution and phased removal of methadone is an acceptable method of detoxification for patients and therapists.

The Methadone assay contained within the Oral Fluid Drug Screen Device yields a positive result when the Methadone concentration in oral fluids exceeds 30 ng/mL.

#### OPIATE (OPI)

The drug class opiates refers to any drug that is derived from the opium poppy, including naturally occurring compounds such as morphine and codeine and semi-synthetic drugs such as heroin. Opiates act to control pain by depressing the central nervous system. The drugs demonstrate addictive properties when used for sustained periods of time; symptoms of withdrawal may include sweating, shaking, nausea and irritability. Opiates can be taken orally or by injection routes including intravenous, intramuscular and subcutaneous; illegal users may also take the intravenously or by nasal inhalation. Using an immunoassay cutoff level of 40 ng/mL, codeine can be detected in the oral fluid within 1 hour following a single oral dose and can remain detectable for 7-21 hours after the dose. 6-monoacetyl-morphine (6-MAM) is found more prevalently in oral fluid and is a metabolic product of heroin. Morphine is the major metabolic product of codeine and heroin and is detectable for 24-48 hours after an opiate dose. The Opiates assay contained within the Oral Fluid Drug Screen Device yields a positive result when the concentration of Morphine in oral fluid exceeds the 40 ng/mL cut-off level.

#### OXYCODONE (OXY)

Oxycodone is a semi-synthetic opioid with a structural similarity to codeine. The drug is manufactured by modifying thebaine, an alkaloid found in opium poppy. Oxycodone, like all opiate agonists, provides pain relief by acting on opioid receptors in the spinal cord, brain, and possibly directly in the affected tissues. Oxycodone is prescribed

for the relief of moderate to high pain under the well-known pharmaceutical trade names of OxyContin®, Tylox®, Percodan® and Percocet®. While Tylox, Percodan and Percocet contain only small doses of oxycodone hydrochloride combined with other analgesics such as acetaminophen or aspirin, OxyContin consists solely of oxycodone hydrochloride in a time-release form.

The Oxycodone assay contained within the Oral Fluid Drug Screen Device yields a positive result when the concentration of oxycodone in oral fluid exceeds 50 ng/mL.

### PHENCYCLIDINE (PCP)

Phencyclidine, the hallucinogen commonly referred to as Angel Dust, can be detected in saliva as a result of the exchange of the drug between the circulatory system and the oral cavity. In a paired serum and saliva sample collection of 100 patients in an Emergency Department, PCP was detected in the saliva of 79 patients at levels as low as 2 ng/mL and as high as 600 ng/mL.

The Phencyclidine assay contained within the Oral Fluid Drug Screen Device yields a positive result when the Phencyclidine concentration in oral fluids exceeds 10 ng/mL.

### MARIJUANA (THC)

Tetrahydrocannabinol, the active ingredient in the marijuana plant (cannabis sativa), is detectable in saliva shortly after use. The detection of the drug is thought to be primarily due to the direct exposure of the drug to the mouth (oral and smoking administrations) and the subsequent sequestering of the drug in the buccal cavity. Historical studies have shown a window of detection for THC in saliva of 3 up to 14 hours after drug use.

The Marijuana assay contained within the Oral Fluid Drug Screen Device yields a positive result when the 11-nor- $\Delta^9$ -THC-9-COOH concentration exceeds 12 ng/mL.

### ALCOHOL (ALCO)

Alcohol intoxication can lead to loss of alertness, coma, death and as well as birth defects.

The B.A.C. at which a person becomes impaired is variable. The United States Department of Transportation (DOT) has established a B.A.C. of 0.02% (20mg/dL) as the cut-off level at which an individual is considered positive for the presence of alcohol.

## PRINCIPLE

### ORAL FLUID DRUG SCREEN DEVICE

The Oral Fluid Drug Screen Device is an immunoassay based on the principle of competitive binding. Drugs that may be present in the oral fluid specimen compete against their respective drug conjugate for binding sites on their specific antibody. During testing, a portion of the oral fluid specimen migrates upward by capillary action. A drug, if present in the oral fluid specimen below its cut-off concentration, will not saturate the binding sites of its specific antibody. The antibody will then react with the drug-protein conjugate and a visible colored line will show up in the test line region of the specific drug strip. The presence of drug above the cut-off concentration in the oral fluid specimen will saturate all the binding sites of the antibody. Therefore, the colored line will not form in the test line region. A drug-positive oral fluid specimen will not generate a colored line in the specific test line region of the strip because of drug competition, while a drug-negative oral fluid specimen will generate a line in the test line region because of the absence of drug competition. To serve as a procedural control, a colored line will always appear at the control line region, indicating that proper volume of specimen has been added and membrane wicking has occurred.

### ALCOHOL TEST

A pad coated with enzymes turns to color shades of green and blue on contact with alcohol in the oral fluids. The alcohol pad employs a solid phase chemistry which uses the following highly specific enzymatic reaction:

During testing, oral fluid is collected on the alcohol pad and saturates the alcohol pad. If no alcohol is present in the oral fluid, the alcohol pad remains colorless (remains white or cream color) because there is no alcohol in the oral fluid to react with enzymes to start the color reaction. If alcohol is present in the oral fluid, the alcohol pad changes to green or blue color because the alcohol reacts with alcohol oxidase to produce aldehyde and peroxide. The peroxide reacts with peroxidase in the presence of hydrogen donor to produce a blue color. Therefore, the presence of green to blue color at the alcohol pad window indicates a presumptive positive result for alcohol.

## WARNINGS AND PRECAUTIONS

1. Do not use after the expiration date.
2. The Oral Fluid Drug Screen Device should remain in the sealed pouch until use.
3. Saliva is not classified as a biological hazard unless derived from a dental procedure.
4. The test device is for single use.
5. Any serious adverse events involving the device should be reported to the manufacturer and the competent authority of the member country where the user and/or patient is located.
6. The used collector and device should be discarded according to state and local regulations.

## CONTENT OF THE PACKAGE

### INCLUDED IN PACKAGE:

- User Instruction
- Multi 11 + ALC Oral Fluid Test (inside foil pouch)
- Desiccants

### NOT INCLUDED IN PACKAGE:

- Timer

## STORAGE CONDITIONS AND VALIDITY

Store as packaged in the sealed pouch at 4-30°C (39-86°F). The test is stable through the expiration date printed on the sealed pouch. The test devices must remain in the sealed pouch until use. **DO NOT FREEZE.** Do not use beyond the expiration date. The test should be used within 30 minutes after being taken out from the aluminum foil bag.


## SPECIMEN REQUIREMENTS

The oral fluid specimen should be collected using the collector provided with the kit. Follow the detailed Directions for Use below. No other collection devices should be used with this assay. Oral fluid collected at any time of the day may be used.

## HOW TO DO THE TEST?

Allow the test device to reach room temperature [15-30°C (59-86°F)] prior to testing. Do not place anything in the mouth including food, drink, gum, or tobacco products for at least 10 minutes prior to collection of oral fluid specimen.

1. Bring the pouch to room temperature before opening it. Remove the test from the sealed pouch and use it as soon as possible.
2. Remove the test device from the sealed pouch and screw the Collector Cap counterclockwise to pull out the entire collection stick with Sponge from the Collection Chamber. (Step 1)
3. Insert the sponge end of the collection stick into the mouth. Close the mouth and gently chew the sponge for saliva excretion. Soak the sponge with saliva in the mouth and swab the inside of the mouth and tongue to collect oral fluid for a total of 3 minutes until the sponge becomes completely soft and fully saturated with saliva. No hard spots should be felt on the sponge when saturated. (Step 2)
4. Remove the sponge from the mouth. With gentle pressure, place the collection stick with the saturated sponge into the Collection Chamber. (Step 3)
5. Screw the Collector Cap clockwise to secure the cap and start the timer. Keep the test tube vertically until you begin to read the test results. (Step 4)
6. Mark the patient information on the test device. Peel off the label to read test results. Wait for the color line(s) to appear on the test strips. (Step 5)
7. **DRUG TEST:** Read results at 10 minutes. **DO NOT** read results after 15 minutes. (Step 6) **ALCOHOL TEST:** Read results at 2 minutes.

## READING THE RESULTS - DRUG TEST

(See image above, Step 6)

### POSITIVE (+)

One color line appears in the control region (C). No line appears in the test region (T).

This positive result indicates that the drug concentration is above the detectable level.

### NEGATIVE (-)

Two lines appear\*. One color line should be in the control region (C), and another apparent color line adjacent should be in the test region (T). This negative result indicates that the drug concentration is below the detectable level.

\*NOTE: The shade of the color in the test line region (T) will vary, but it should be considered negative whenever there is even a faint distinguishable color line.

### INVALID


Control line fails to appear in 10 minutes. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test using a new test device. If the problem persists, discontinue using the lot immediately and contact your supplier.

**NOTE:** Each test strip must be looked at individually. Each line may vary in color and darkness. **(DO NOT)** compare lines within the same test strip or between different test strips).

A positive test result does not always mean a person took illegal drugs and a negative test result does not always mean a person did not take illegal drugs. There are several factors that influence the reliability of drug tests. Certain drugs of abuse tests are more accurate than others.

**IMPORTANT:** The result you obtained is called preliminary for a reason. The sample must be tested by laboratory to determine if a drug of abuse is actually present.

## READING THE RESULTS - ALCOHOL TEST

### ALCOHOL NEGATIVE RESULT

The alcohol pad shows no color change (remains white or cream colored); it should be interpreted as a negative result (no alcohol present). A result where the outer edges of the alcohol pad produces a slight color but the majority of the pad remains colorless should be repeated to ensure complete saturation of the alcohol pad with oral fluid. If the second result is the same, the results should be interpreted as being negative (no alcohol present).

### ALCOHOL PRESUMPTIVE POSITIVE RESULT

The alcohol test produces a color change to green to blue in the presence of salivary alcohol 0.02% B.A.C. or higher. At higher alcohol concentration near 0.30% B.A.C., the color may change to a dark blue-gray.

### WHAT IS A FALSE POSITIVE TEST?

The definition of a false positive test would be an instance where the test result from the LumxedIQ One-Step Drug of Abuse Test is positive, even though the initial target drug is not present in the sample. The most common causes of a false positive test are cross reactants. Certain foods and medicines, dietary plan drugs and nutritional supplements may also cause a false positive test result with this product.

### WHAT IS A FALSE NEGATIVE TEST?

The definition of a false negative test is that the initial target drug is present in the oral fluid sample but is not detected by the LumxedIQ One-Step Drug of Abuse Test. A diluted sample or a sample contaminated with a substance could cause a false negative result.

### QUALITY CONTROL

A procedural control is included in the test. A color line appearing in the control region (C) is considered an internal procedural control. It confirms sufficient specimen volume, adequate membrane wicking and correct procedural technique.

### TEST LIMITATIONS

- The Oral Fluid Drug Screen Device provides only a qualitative, preliminary analytical result. A secondary, analytical method must be used to obtain a confirmed result. Gas chromatography/mass spectrometry (GC/MS) or Liquid chromatography/tandem mass spectrometry (LC/MS/MS) are preferred confirmatory methods.
- A positive test result does not indicate the concentration of drug in the specimen or the route of administration.
- A negative result may not necessarily indicate a drug-free specimen. Drug may be present in the specimen below the cutoff level of the assay.
- The test has been developed for testing saliva sample only. No other fluids have been evaluated. DO NOT use this device to test anything but saliva.

## PERFORMANCE CHARACTERISTICS

### ANALYTICAL SENSITIVITY

A negative saliva pool was spiked with drugs to target concentrations of  $\pm$  50% cut-off and  $\pm$  25% cut-off and tested with the Oral Fluid Drug Screen Device. The results are summarized below.

Drug concentration Cut-off Range	n	AMP50		BARB50		BZO50		COC20	
		-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	50	50/0
-50% Cut-off	30	30	0	30	0	30	0	50	50/0
-25% Cut-off	30	28	2	28	2	28	2	50	50/0
Cut-off	30	13	17	13	17	13	17	50	50/0
25% Cut-off	30	4	26	4	26	4	26	50	2/48
50% Cut-off	30	0	30	0	30	0	30	50	0/50

Drug concentration Cut-off Range	n	mAMP50		MDMA50		MTD30		OPI40	
		-	+	-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0	30	0
-25% Cut-off	30	29	1	29	1	29	1	27	3
Cut-off	30	16	14	5	25	16	14	18	12
25% Cut-off	30	7	23	0	30	7	23	3	27
50% Cut-off	30	0	30	0	30	0	30	0	30

Drug concentration Cut-off Range	n	OXY50		PCP10		THC12	
		-	+	-	+	-	+
0% Cut-off	30	30	0	30	0	30	0
-50% Cut-off	30	30	0	30	0	30	0
-25% Cut-off	30	28	2	28	2	29	1
Cut-off	30	12	18	13	17	12	18
25% Cut-off	30	3	27	4	26	3	27
50% Cut-off	30	0	30	0	30	0	30

For the alcohol test saliva was obtained by rinsing with positive ethanol control solutions at various B.A.C. (0.02%, 0.08%, 0.15%, 0.30%). Negative saliva was used to test at 0% concentration. For each concentration, a total of 15 tests were performed to validate the test performance. The results of the Oral Fluid Drug Screen Device are summarized below:

Test	Total # of Test / Concentration	B.A.C.									
		0.00%		0.02%		0.08		0.15%		0.30%	
		-	+	-	+	-	+	-	+	-	+
Alcohol	15	15	0	1	14	0	15	0	15	0	15

## CLINICAL PERFORMANCE

### AMPHETAMINE (AMP)

AMP (Amphetamine 50ng/ml) (LC/MS/MS values calibrated to D-Amphetamine): In this study, one hundred and five negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	21	31
Negative	23	8	20	1	0

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.11%	90.06%	99.67%
Diagnostic specificity	98.08%	89.88%	99.66%
Total coincidence rate	98.10%	93.32%	99.48%

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.11%	90.06%	99.67%
Diagnostic specificity	98.08%	89.88%	99.66%
Total coincidence rate	98.10%	93.32%	99.48%

### BARBITURATES (BARB)

BARB (Barbiturates 50ng/ml) (LC/MS/MS values calibrated to Butalbital): In this study, one hundred and three negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	19	32
Negative	22	9	20	1	0

### METHAMPHETAMINES (MAMP)

mAMP (Methamphetamines 50ng/ml) (LC/MS/MS values calibrated to D-Methamphetamine): In this study, one hundred and four negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	19	32
Negative	22	9	20	1	0

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.08%	89.88%	99.66%
Diagnostic specificity	98.04%	89.70%	99.65%
Total coincidence rate	98.06%	93.19%	99.47%

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.08%	89.88%	99.66%
Diagnostic specificity	98.08%	89.88%	99.66%
Total coincidence rate	98.08%	93.26%	99.47%

### BENZODIAZEPINES (BZO)

BZO (Benzodiazepines 50ng/ml) (LC/MS/MS values calibrated to Oxazepam): In this study, one hundred and five negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	21	31
Negative	23	8	20	1	0

### METHYLENEDIOXYMETHAMPHETAMINE (MDMA)

MDMA (Methylenedioxyamphetamine 50ng/ml) (LC/MS/MS values calibrated to (±)-3,4-Methylenedioxyamphetamine): In this study, one hundred and four negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	21	30
Negative	22	9	20	1	0

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.08%	89.88%	99.66%
Diagnostic specificity	98.08%	89.88%	99.66%
Total coincidence rate	98.08%	93.26%	99.47%

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.08%	89.88%	99.66%
Diagnostic specificity	98.08%	89.88%	99.66%
Total coincidence rate	98.08%	93.26%	99.47%

### COCAINE (COC)

COC (Cocaine 20ng/ml) (LC/MS/MS values calibrated to Benzoylcegonine): In this study, one hundred and two negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	20	31
Negative	20	9	20	1	0

### METHADONE (MTD)

MTD (Methadone 30ng/ml) (LC/MS/MS values calibrated to Methadone): In this study, one hundred and four negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	19	32
Negative	22	9	20	1	0

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.08%	89.88%	99.66%
Diagnostic specificity	98.08%	89.88%	99.66%
Total coincidence rate	98.08%	93.26%	99.47%

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.08%	89.88%	99.66%
Diagnostic specificity	98.08%	89.88%	99.66%
Total coincidence rate	98.08%	93.26%	99.47%

### OPIATES (OPI)

OPI (Opiates 40ng/ml) (LC/MS/MS values calibrated to Morphine): In this study, one hundred and six negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	2	20	32
Negative	22	9	20	1	0

### MARIJUANA (THC)

THC (Marijuana 12ng/ml) (LC/MS/MS values calibrated to 11-nor- $\Delta^8$ -THC-9-COOH): In this study, one hundred and six negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	2	20	32
Negative	22	9	20	1	0

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.11%	90.06%	99.67%
Diagnostic specificity	96.23%	87.25%	98.96%
Total coincidence rate	97.17%	92.01%	99.03%

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.11%	90.06%	99.67%
Diagnostic specificity	96.23%	87.25%	98.96%
Total coincidence rate	97.17%	92.01%	99.03%

### OXYCODONE (OXY)

OXY (Oxycodone 50ng/ml) (LC/MS/MS values calibrated to Oxycodone): In this study, one hundred and two negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	19	31
Negative	21	9	20	1	0

### ALCOHOL (ALCO)

Alcohol (Alcohol 20mg/dL) (GC values calibrated to alcohol): In this study, one hundred and seven negative and positive saliva samples were tested and compared with GC. The results are summarized below:

WHPM Result	Reference Method: Concentration $\geq$ 20mg/dL	Reference Method: Concentration < 20mg/dL	Total
Positive	53	1	54
Negative	1	52	53
Total	54	53	107

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.04%	89.70%	99.65%
Diagnostic specificity	98.04%	89.70%	99.65%
Total coincidence rate	98.04%	93.13%	99.46%

/	Coincidence rate	Lower limit of 95% confidence interval	Upper limit of 95% confidence interval
Diagnostic sensitivity	98.15%	90.23%	99.67%
Diagnostic specificity	98.11%	90.06%	99.67%
Total coincidence rate	98.13%	93.44%	99.49%

### PHENCYCLIDINE (PCP)

PCP (Phencyclidine 10ng/ml) (LC/MS/MS values calibrated to Phencyclidine): In this study, one hundred and six negative and positive saliva samples were tested and compared with LC/MS/MS. The results are summarized below:

WHPM Result	Drug-Free	Less Than Half the Cutoff Concentration by GC/MS Analysis	Near Cutoff Negative (Between 50% below the Cutoff and the Cutoff Concentration)	Near Cutoff Positive (Between the Cutoff and 50% above the Cutoff Concentration)	High Positive (Greater Than 50% above the Cutoff Concentration)
Positive	0	0	1	21	30
Negative	21	9	21	1	0

### CROSS REACTIVITY

The following tables list the concentration of compounds above which the Oral Fluid Drug Screen Device identified positive results at a read time of 10 minutes.

### AMPHETAMINE (AMP)

Drug	Concentration (ng/mL)
D-Amphetamine	50
DL-Amphetamine	125
$\beta$ -Phenylethylamine	4,000
( $\pm$ )-3,4-Methylenedioxyamphetamine (MDA)	150
L-Amphetamine	4,000
p-Hydroxyamphetamine	800
Tryptamine	1,500
Tyramine	1,000

## BARBITURATES (BARB)

Drug	Concentration (ng/mL)
Alphenol	300
Amobarbital	150
Aprobarbital	300
Butabarbital	200
Butethal	75
Cyclopentobarbital	100
Pentobarbital	600
Phenobarbital	300
Secobarbital	100

## BENZODIAZEPINES (BZO)

Drug	Concentration (ng/mL)
a-Hydroxylprazolam	1,260
Alprazolam	40
Bromazepam	400
Chlordiazepoxide	780
Chlordiazepoxide HCl	390
Clobazam	100
Clonazepam	785
Clorazepate Dipotassium	195
Delorazepam	1,560
Desalkylflurazepam	390
Diazepam	195
Estazolam	2,500
Flunitrazepam	385
(±) Lorazepam	1,560
RS-Lorazepam glucuronide	160
Midazolam	12,500
Nitrazepam	95
Norchlordiazepoxide	200
Nordiazepam	390
Oxazepam	50
Temazepam	20
Triazolam	2,500

## COCAINE (COC)

Drug	Concentration (ng/mL)
Benzoyllecgonine	20
Cocaine HCl	20
Cocaine ethylene	25
Ecgonine HCl	1,500
Ecgonine methyl ester	12,500

## METHAMPHETAMINES (MAMP)

Drug	Concentration (ng/mL)
D-Methamphetamine	50
(1R,2S) - (-) Ephedrine	400
Fenfluramine	60,000
Methoxyphenamine	25,000
(±)-3,4-Methylenedioxyamphetamine (MDMA)	50
p-Hydroxymethamphetamine	400
L-Phenylephrine	4,000
Procaine	2,000

## METHYLENEDIOXYMETHAMPHETAMINE (MDMA)

Drug	Concentration (ng/mL)
(±)-3,4-Methylenedioxyamphetamine	50
Dobutamine Hydrochloride Injection	60,000
p-Hydroxymethamphetamine	15,000
(±)-3,4-Methylenedioxyamphetamine	1,500

## METHADONE (MTD)

Drug	Concentration (ng/mL)
Methadone	30
Doxylamine	12,500

## OPIATES (OPI)

Drug	Concentration (ng/mL)
Morphine	40
Bilirubin	3,500
Codeine	10
Diacetylmorphine (Heroin)	50
Ethylmorphine	24
Hydrocodone	100
Hydromorphone	100
Levorphanol	400
6-Monoacetylmorphine	25
Morphine 3-β-D-Glucuronide	50
Nalorphine	10,000
Normorphine	12,500
Norcodeine	1,500
Oxycodone	25,000
Oxymorphone	25,000
Thebaine	1,500

## OXYCODONE (OXY)

Drug	Concentration (ng/mL)
Oxycodone	50
Codeine	25,000
Dihydrocodeine	6,250
Ethylmorphine	12,500
Hydrocodone	1,000
Hydromorphone	6,250
Oxymorphone	1,000
Thebaine	25,000

## PHENCYCLIDINE (PCP)

Drug	Concentration (ng/mL)
Phencyclidine	10
Tetrahydrozoline	50,000

## MARIJUANA (THC)

Drug	Concentration (ng/mL)
11-nor-Δ <sup>9</sup> -THC-9-COOH	12
Cannabinol	3,000
Δ <sup>9</sup> -THC	75
Δ <sup>8</sup> -THC	15

## ALCOHOL TEST

The alcohol test will react with methyl, ethyl, and allyl alcohols, but it will not react with alcohols having 5 or more carbons like glycine, glycerol, and serine. This property is a result of specificity of the alcohol oxidase enzyme extracted from yeast.

## INTERFERING SUBSTRATE

A study was conducted to determine the cross-reactivity of the test with compounds spiked into drug-free PBS stock. The following compounds demonstrated no false positive results on the Oral Fluid Drug Screen Device when tested with concentrations up to 100 µg/mL.

Acebutolol Hydrochloride	Isoproterenol Hydrochloride Injection
Acepromazine-d6 hydrochloride	Ketoprofen
Acetaminophen	Labetalol Hydrochloride
N-Acetylprocainamide	Loperamide Hydrochloride
Acetophenetidin	Loxapine succinate salt
Amoxicillin	Maprotiline hydrochloride
Ampicillin	Meperidine
Amitriptyline Hydrochloride Tablets	Meprobamate
R-(-)-Apomorphine Hydrochloride	Methylphenidate
Aspirin	Methprylon
Aspartame	Morphine-3β-D-glucuronide
L-Ascorbic Acid	Nalidixic acid
Atropine	Naproxen
Benzphetamine HCL	Naloxone
Benzilic acid	Naltrexone hydrochloride
Benzoylcegonine	Nicotinamide (vitamin B3)
Benzoic Acid	Nimesulide
Bilirubin, MixedIsomers	Nifedipine
Brompheniramine maleate	Nordoxepin hydrochloride
Buspirone hydrochloride	Norethisterone Tablets
Cannabidiol	Noscapine HCL Hydrate
Chlorothiazide	Noroxymorphone HCL
Chloroquine diphosphate	Nylidrin hydrochloride
Chlorpheniramine Maleate	(±)-Octopamine HCL
Chlorpromazine HCL	Oxalic Acid
Chloramphenicol	Oxolinic Acid
Chloral Hydrate	Oxymetazoline Hydrochloride
Cholesterol	Papaverine hydrochloride
Chlorothiazide	Pentazocine
Clomipramine Hydrochloride Tablets	Perphenazine
Clonidine hydrochloride	Penicillin G Sodium salt
Cortisone	Phenelzine sulfate salt
Creatinine	Phentermine HCL
Dextromethorphan	Phenylpropanolamine hydrochloride
Diclofenac Sodium	Prednisolone
Dicyclomine	Prednisone Acetate Tablets
Diffunisal	Promazine hydrochloride
Digoxin	Promethazine
4-Dimethylaminoantipyrine	Propranolol Hydrochloride
5,5-Diphenylhydantoin	Pseudoephedrine HCL
Diphenhydramine	Quinine
Dopamine Hydrochloride Injection	Quinidine
Ecgonine methylester	Quinacrine dihydrochloride
Efavirenz	Ranitidine Hydrochloride Tablets
Emetine dihydrochloride hydrate	Salicylic Acid
(-)-Epinephrine	Serotonin HCl
Erythromycin	Sertraline HCl
Estradiol	Sulfamethazine
Estrone-3-sulfate potassium salt	Sulindac
Ethyl-p-aminobenzoate	Tetracycline
Fenopropfen calcium salt hydrate	Tetrahydrocortisone 3-(β-D-glucuronide)
Furosemide	Theophylline (Sustained-release Tablets)
Gentisic acid	Thioridazine
D-Glucuronic acid	Thiamine, (Vitamin B1)Tablets)HCL

Glutethimide	L-Thyroxine
Guaifenesin Guaiaacol glyceryl ether	Tolbutamide
Hemoglobin porcine	Trimethoprim
Hippuric Acid	Trazodone hydrochloride
Hydralazine hydrochloride	Triamterene
a-Hydroxyhippuric acid	Trimipramine
21-Hydroxy progesterone	Tryptamine
Hydrocortisone	DL-Tryptophan
Hydrochlorothiazide	Trans-2-phenylcyclopropylamine hydrochloride
(+/-)-4-Hydroxyamphetamine HCL	D/L-Tyrosine
Ibuprofen	Tyramine
Imipramine HCL	Uric Acid
Iprazid	Verapamil Hydrochloride
Isosuprine hydrochloride	Zomepirac sodium salt

## ALCOHOL TEST

The following substances may interfere with the Oral Fluid Drug Screen Device when using samples other than oral fluid:


1. Agents which enhance color development: Peroxides and strong oxidizers.
2. Agents which inhibit color development.
3. Reducing agents: such as ascorbic acid, tannic acid, pyrogallol, mercaptanals and tosylates, oxalic acid, uric acid, bilirubin, L-dopa, L-methyldopa, and methampyrone, etc.

The above-named substances do not normally appear in sufficient quantity in oral fluid to interfere with the test. However, care must be taken that they are not introduced into the mouth during the 10 minutes period preceding the test.

## BIBLIOGRAPHY OF SUGGESTED READING

1. Moolchan, E., et al, "Saliva and Plasma Testing for Drugs of Abuse: Comparison of the Disposition and Pharmacological Effects of Cocaine", Addiction Research Center, IRP, NIDA, NIH, Baltimore, MD. As presented at the SOFT-TIAFT meeting October 1998.
2. Kim, I, et al, "Plasma and oral fluid pharmacokinetics and pharmacodynamics after oral codeine administration", Clin Chem, 2002 Sept.; 48 (9), pp 1486-96.
3. Schramm, W. et al, "Drugs of Abuse in Saliva: A Review," J Anal Tox, 1992 Jan-Feb; 16 (1), pp 1-9.
4. McCarron, MM, et al, "Detection of Phencyclidine Usage by Radioimmunoassay of Saliva," J Anal Tox. 1984 Sep-Oct; 8 (5), pp 197-201.
5. Tietz NW. Textbook of Clinical Chemistry. W.B. Saunders Company. 1986; 1735.

## SYMBOLS

	Do not reuse		Catalog number
	Temperature limitations		Consult instructions before use
	Manufacturer		Authorized Representative in the European Union
	Batch code		Use by
	Contains sufficient for <n> tests		

**DISTRIBUTED BY:**

LuxmediQ Deutschland UG  
Alonaer Strasse 27, 10555 Berlin, Germany

W.H.P.M. Bioresearch & Technology Co., Ltd.  
No. 2 Zhongxin Street, Louzizhuang, Jinzhan Xiang, Chaoyang District, Beijing,  
10018 P.R. China  
Tel: (+86)-10-84391888  
Fax: (+86)-10-84391888-117  
www.whpm.com.cn

Effective date: 2021-12-10  
Version 1.0

**FRAGEN | QUESTIONS?**

E-mail: [service@luxmediq.de](mailto:service@luxmediq.de)

Tel: +49 6401 9699 390

Web: [www.luxmediq.de](http://www.luxmediq.de)