LuxmedIQ

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- Stewart D J, Inaba T, Lucasen M, et al. Cocaine metabolism: cocaine and norcocaine hydrolysis by liver and serum esterases [J]. *Clinical Pharmacology & Therapeutics*, 1979, 25(4):464.
- Ambre J. The urinary excretion of cocaine and metabolites in humans: a kinetic analysis of published data[J]. *Journal of Analytical Toxicology*, 1985, 9(6):241.
- Hawks R L, Chiang C N. Urine testing for drugs of abuse[J]. *Nida Research Monograph*, 1986, 73(2):1.

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LuxmedIQ

LuxmedIQ Single/Multi Rapid Drug Test (Urine)

The LuxmedIQ One-Step Drug of Abuse Test is a rapid test for the qualitative detection of drugs or their breakdown products in human urine. **Please read these instructions for use carefully before testing, as they contain important information about what you need to consider when testing.**


NOTE: The test may yield preliminary positive results when prescription drugs are ingested at prescribed doses. It is not intended to distinguish between prescription use and abuse of any drug. There are no uniformly recognized cutoff concentration levels for any drug in urine. The test provides only preliminary test results. A more specific alternative chemical method must be used to obtain a confirmed analytical result. Gas chromatography/mass spectrometry (GC/MS) is the preferred confirmatory method. Clinical consideration and professional judgment should be exercised with any drug of abuse test result, particularly when the preliminary result is positive. This is not a screening device to monitor prescription medication. This test does not provide any basis for medical interpretation.

WHAT IS A ONE-STEP DRUG OF ABUSE TEST?

After drug consumption, the drug is metabolized by the body and excreted unchanged, or as a breakdown product, via the urine. Due to this metabolic process, the concentration of the drug and/or its degradation products in the urine initially increases, but levels off over time and then decreases again. The LuxmedIQ One Step Drug Test is a rapid test for the qualitative detection of drugs or their breakdown products in human urine.

WHAT IS THE CUTOFF AND THE DETECTION WINDOW?

The cutoff concentration of a drug test is the concentration of a substance in the urine above which the test provides a positive result. The drug test therefore provides a preliminary positive result if the concentration of the drug and/or its breakdown products in the urine reaches or exceeds the cut-off concentration of the test. If the concentration of the drug or its degradation products in the urine is below the cut-off concentration, the test provides a negative result. The following image illustrates how the test works.

The following table provides an overview of the earliest time a drug can be detected in urine after consumption (Minimum Detection Time) and how long it can be detected in urine (Maximum Detection Time). The time in which the drug can be detected with the test is called the detection window.

NOTE: The following table contains non-binding reference values. The exact detection window depends on various individual factors (amount and frequency of drug consumed, individual metabolism, the drug itself, etc.).


Drug	Minimal Detection Time	Maximal Detection Time
Alcohol (ALC)	1 hour after consumption	1-3 days
Amphetamines (AMP)	4-6 hours after consumption	1-3 days (highly dependent on the pH value of the urine)
Barbiturates (BAR)	4-6 hours after consumption	1 day (short-acting) to 20 days (long-acting)

Buprenorphine (BUP)	4-6 hours after consumption	2-5 days
Benzodiazepines (BZO)	4-6 Hours after consumption	2-3 days, several weeks in the case of chronic consumption
Cocaine (COC)	2-6 Hours after consumption	2-4 days
Cotinine (COT)	4 Stunden circa	2-3 days
EDDP (methadone degradation product)	4-6 Hours after consumption	2-7 days
Ethyl glucuronide (ETG)	1-2 Hours after consumption	5 days, with heavy consumption
Fentanyl (FEN)	1-2 Hours after consumption	1-2 days
Ketamine (KET)	2-6 Hours after consumption	2-4 days
Ecstasy (MDMA)	4-6 Hours after consumption	1-3 days (highly dependent on the pH value of the urine)
Methamphetamine (MET, mAMP)	4-6 Hours after consumption	1-3 days
Opiate-Morphines (MOP/OPI)	2-6 Hours after consumption	2-4 days
Methylphenidate (MPD)	4-6 Hours after consumption	1-2 days
Metadone (MTD)	4-6 Hours after consumption	2-5 days
Oxycodone (OXY)	4-6 Hours after consumption	1-3 days
Phencyclidine (PCP)	4-6 Hours after consumption	2-4 days
Propoxyphene (PPX)	4-6 Hours after consumption	1-2 days
K2 (synth. cannabinoids)	1-3 Stunden nach dem Konsum	1-3 days
Antidepressants (TCA)	4-6 Hours after consumption	2-3 days
Cannabinoids (THC)	1-3 Hours after consumption	2-3 days, with regular consumption several weeks
Tramadol (TRA)	4-6 hours after consumption	2-6 days

PRINCIPLE

The LuxmedIQ One-Step Drug of Abuse Test is an immunoassay. During testing, the urine specimen migrates upward on the test strip. A drug-positive urine specimen will not generate a colored line in the specific test line region of the strip, while a drug-negative urine specimen will generate a line in the test line region. A colored line will always appear in the control line region, indicating that the proper volume of specimen has been added.

The following image shows the structure of the multi-test (cassette) on the left and the single test (test strips) on the right. The multi-test consists of several individual tests (test strips), which are enclosed in a test cassette. The functionality of both test types is therefore identical.

WARNINGS AND PRECAUTIONS

- For external (outside of the human body) use only.
- For single use. Discard after first use.
- Do not use the test if the pouch is punctured or not well sealed.
- Do not use after the expiration date.
- Keep out of the reach of children.
- The used test cup and urine specimen should be discarded according to federal, state and local regulations.

- The used test device and urine specimen should be discarded according to federal, state and local regulations.

PACKAGE CONTENTS

INCLUDED IN THE PACKAGE:

- User instructions
- Test strip (individually sealed in a foil pouch)
- Desiccants

NOT INCLUDED IN THE PACKAGE:

- Watch, timer or clock
- Disposable gloves
- Collection cup for the urine sample

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 36°F - 86°F (2°C - 30°C). The test is stable through the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. If the test remains unused for more than 10 minutes while removed from the sealed pouch, results may be erroneous. Keep away from direct sunlight, moisture, and heat. **DO NOT FREEZE.** Do not use beyond the expiration date.

WHEN TO COLLECT THE URINE SAMPLE FOR THE TEST?

Urine collected at any time of the day can be used.


TEST PREPARATION

- When you are ready to test, remove the strip from the sealed foil pouch.
- Notice the colored tape on the strip indicates the name of the drug you are testing for.
- Fill a sample cup with a fresh urine sample. **DO NOT FREEZE.**

HOW TO DO THE TEST?

- Do not take the test strip / multitest (cassette) out of the sealed pouch until immediately before performing the test. Remove the protective cap from the multitest.
- Immerse the test strip / multitest in the urine sample with the urine absorption area pointing downwards. Do not immerse the test deeper than the MAX (maximum) line. When using the multitest, do not allow any urine to touch the plastic housing of the cassette. Remove the test strip / multitest from the urine sample after at least 15 seconds. Lay the test flat on a non-absorbent, clean surface. Replace the protective cap on the Multitest.
- Wait for 5 minutes (start the timer immediately after the test has been removed from the urine sample) and read the result. **DO NOT** read the result after more than 10 minutes.

NOTE: Results read after over 10 minutes may not be accurate and should not be read.

READING THE RESULTS

PRELIMINARY POSITIVE (+)

If a line appears in the control area (C) but NO line appears in the drug test area (T), then this indicates a preliminary positive result for the corresponding drug.

NEGATIVE (-)

If a line appears in both the control area (C) and the drug test area (T), then this indicates a negative result for the corresponding drug, regardless of how dark or light the line may appear.

INVALID

If at 5 minutes, NO line appears in the control area (T), then the result is invalid.

NOTE: Each test strip must be looked at individually. Each line may vary in color and darkness. (DO NOT compare lines within the same test strip or between different test strips).

A positive test result does not always mean a person took illegal drugs and a negative test result does not always mean a person did not take illegal drugs. There are several factors that influence the reliability of drug tests. Certain drugs of abuse tests are more accurate than others.

IMPORTANT: The result you obtained is called preliminary for a reason. The sample must be tested by laboratory to determine if a drug of abuse is actually present.

WHAT IS A FALSE POSITIVE TEST?

The definition of a false positive test would be an instance where the test result from the LuxmedIQ One-Step Drug of Abuse Test is positive, even though the initial target drug is not present in the sample. The most common causes of a false positive test are cross reactants. Certain foods and medicines, dietary plan drugs and nutritional supplements may also cause a false positive test result with this product.

WHAT IS A FALSE NEGATIVE TEST?

The definition of a false negative test is that the initial target drug is present in the urine sample but is not detected by the LuxmedIQ One-Step Drug of Abuse Test. If the sample is diluted, or if the sample is contaminated with a substance this could cause false negative results.

BIBLIOGRAPHY OF SUGGESTED READINGS

- Stewart D J, Inaba T, Lucassen M, et al. Cocaine metabolism: cocaine and norcocaine hydrolysis by liver and serum esterases[J]. *Clinical Pharmacology & Therapeutics*, 1979, 25(4):464.
- Amre J. The urinary excretion of cocaine and metabolites in humans: a kinetic analysis of published data[J]. *Journal of Analytical Toxicology*, 1985, 9(6):241.
- Hawks R L, Chiang C N. Urine testing for drugs of abuse[J]. *Nida Research Monograph*, 1986, 73(2):1.

SYMBOLS

	Do not reuse		Read instructions before use
	Temperature limitations		Manufacturer
	Batch code		Use by
	Contains sufficient for <n> tests		

LuxmedIQ Deutschland UG
 Altonaer Strasse 27
 10555 Berlin
 Germany

Effective date: 2023-12-06

FRAGEN | QUESTIONS?

E-mail: service@luxmediq.de

Tel: +49 6401 9699 390

Web: www.luxmediq.de