Sofort-Kältekompresse

1. PRODUKTNAME

Kältekompressen

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17. HERSTELLER

JIANGSU INTCO MEDICAL PRODUCTS CO., LTD.
No. 298 Yandunshan Road, Zhenjiang New District, 212132 Jiangsu, P.R. China

18. EU-BEVOLLMÄCHTIGTER

MedNet EC-REP C Ilb GmbH
Borkstrasse 10, 48163 Münster, Germany

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2026-06-23

20. REVISION:

A/2

21. KONTAKTINFORMATIONEN

JIANGSU INTCO MEDICAL PRODUCTS CO., LTD.
No. 298 Yandunshan Road, Zhenjiang New District, 212132 Jiangsu, P.R. China
E-Mail: sunguihua@intco.com

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C128150NA / Contains: urea, water


3. METHOD OF USE

1. Shake the pouch so that the contents collect at the bottom. Locate the inner pouch. Squeeze firmly to rupture the inner fluid pouch.

2. Shake briefly to mix the contents.

3. Wrap the pack in a towel or cloth before use.

4. Apply to the affected area. The recommended application time is 15–20 minutes.

4. INTENDED PURPOSE

Provides cold therapy for the temporary relief of pain, reduction of swelling, and cooling of the affected area. Suitable for use in healthcare facilities or at home. Refer to the Instructions for Use for warnings, precautions, and limitations of use.

5. INTENDED PATIENT POPULATION

The cold pack is suitable for use with most patients. However, extra care should be taken when using it with elderly patients or children. Monitor the skin regularly throughout use. Do not use the cold pack on the following patients:

- Patients with impaired sensation or patients who are unconscious.
- Pregnant women and infants.
- Patients with any of the conditions mentioned in the contraindications below.

6. INTENDED USER

Cold packs must only be applied by an adult who has read and understood the Instructions for Use and is able to follow them.

7. INDICATIONS

- Relieves local pain and swelling caused by sprains, strains, impacts, bruises, and other injuries when used within 48–72 hours after the injury occurs;
- Relieves local pain and discomfort caused by toothache and mosquito bites;
- Provides physical cooling for patients with a high fever;
- Relieves local pain and swelling caused by mild burns;
- Reduces swelling, pain, and bleeding around the wound after surgery.

8. CONTRAINDICATIONS

- Raynaud's phenomenon;
- Cold hypersensitivity;
- Cardiac conditions;
- High blood pressure;
- Do not apply over areas where nerves lie superficially;
- Mid-to-late stage of the wound.

9. UNDESIRABLE SIDE EFFECTS

Avoid pressing the cold pack directly against the skin, as this may cause a cold injury. In most cases, this symptom is transient and typically resolves within a few hours. If the symptom is not alleviated within one day, please seek medical attention.

10. INTENDED CLINICAL BENEFITS

Relieves pain, controls inflammation and swelling, promotes hemostasis, provides physical cooling, and controls postoperative swelling.

11. WARNINGS

- For external use only.
- Do not use the product if the packaging is damaged or leaking, or after the expiry date.
- Do not use on infants.
- Harmful if swallowed. If accidentally swallowed, contact a poison control center or physician immediately.
- If the cold pack breaks and the fluid comes into contact with the skin or eyes, flush with water immediately and call a doctor.
- Individuals with circulatory problems should consult a physician before using this product.
- Pay particular attention when using the product on elderly patients, children, and people with communication barriers.

12. CAUTIONS

- Read the instructions carefully before use.
- Prior to use, store at room temperature.
- Do not apply directly to the skin.
- Always wrap the pack in a towel or cloth before applying it to the skin.
- If the cold pack feels too cold, or if there are noticeable changes in skin color, add additional layers of towel or cloth or remove the cold pack.
- Do not lie on the cold pack or place your full body weight on it.
- Do not open the cold pack.
- Do not reuse by freezing or refrigerating.
- Keep out of reach of children.

13. DISPOSAL

Dispose of the pack in accordance with applicable local, regional, national, or international regulations.

14. STORAGE/TRANSPORTATION CONDITIONS

Store and transport at temperatures between 4 °C and 35 °C, away from direct sunlight.

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17. MANUFACTURER

 JIANGSU INTCO MEDICAL PRODUCTS CO., LTD.
No. 298 Yandunshan Road, Zhenjiang New District, 212132 Jiangsu, P.R. China

18. EU AUTHORIZED REPRESENTATIVE

 MedNet EC-REP C IIb GmbH
Borkstrasse 10, 48163 Münster, Germany

19. EFFECTIVE DATE:

2026-06-23

20. REVISION:

A/2

21. CONTACT INFORMATION

JIANGSU INTCO MEDICAL PRODUCTS CO., LTD.
No. 298 Yandunshan Road, Zhenjiang New District, 212132 Jiangsu, P.R. China
Email: sunguihua@intco.com

For patients, users and third parties in the European Union and in countries with an equivalent regulatory framework in the European Union and in countries with an identical regulatory regime (Regulation (EU) 2017/745 on medical devices): If a serious incident occurs during the use of this device or as a result of its use, please report it to the manufacturer and/or its authorized representative and to your national competent authority.

Contact details for national competent authorities (Vigilance Contact Points) and further information can be found on the following European Commission website:

https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

22. E-IFU:

 <https://www.intcomedical.com/download.html>
Website: www.intcomedical.com

23. SYMBOLS

	Catalog number		Temperature limit
	Consult electronic instructions for use		Batch code
	Caution		Use-by date
	Manufacturer		Date of Manufacture Country of Manufacture
	Importer		Unique device identifier
	Do not use if package is damaged		Do not reuse
	Medical device		CE marking and Notified Body code
	EU authorized representative		Model number
	Keep away from sunlight		Keep dry

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